

## Quickfix detectable - Ref. 7304



### Identification of the substance/preparation and of the company/undertaking

- Product: QuickFix Detectable.
- Company: Plum A/S.

### Composition

- Blue PEplasters.
- Wound pad with detectable aluminium foil placed centrally.
- Individually packed.
- Sterilized with gamma rays.
- Glue: Acryl polymer.
- Water-resistant: Yes.

### First-aid measures

- General: None
- Inhalation: -
- Contact with skin: -
- Contact with Eyes: -
- Ingestion: If necessary, seek medical advice and bring data sheet and label.
- Hazards Identification: None.

### Fire-fighting measures

- Fire extinguisher: All extinguishers can be used.

### Accidental release measures

- Ecological precautions: None.

- Personal precautions: None.

### Handling and storage

- Handling: No special precautions.

- Storage: Best at 0-30°C.

### Exposure controls / personal protection

- None is recommended.

### Physical properties

- Appearance: Blue.
- Size: L 7,2 x B 2,5 cm.
- Form: PE-plaster, water-resistant, detectable.

### Stability and reactivity

- Special properties or risks: No special precautions.

### Toxicological information

- Acute toxicity: -
- Hazards to health: -
- Inhalation: -
- Contact with skin: The product is tested against skin irritation ; Hypoallergenic.
- Contact with eyes: -
- Ingestion: -
- Sensitisation: All contents have been chosen with the purpose of reaching the lowest possible risk for sensitisation.

### Ecological information

- None.

### Disposal

- Remains of the product and the packaging are combustible and do not give off any poisonous matters to the atmosphere.

### Transport information

- The product is no dangerous goods. No special information required regarding transportation.

### Regulatory information

- None.

### Other information

- Application: For use in food industries and other places where detect ability and visibility is important.
- Directions for use: Pull the plaster out of the refill and on the wound.
- Shelf life: min. 30 months.